
Mask

Full Face Mask (F3)

User Manual

CE 0123

Components of Full Face Mask (F3)

A. Frame	H. Oxygen / Pressure Port Connector
B. Forehead pad	I. Anti-asphyxiation valve
C. Cushion	J. Exhaust holes
D. Headgear	K. Elbow "V" (transparent)
E. Connector	L. Elbow "NV1" (milky)
F. Stoppie	M. Elbow "NV2" (blue)
G. Swivel	N. Elbow "NV3" (green)

Full Face Mask (F3)

Thank you for choosing BMC's Full Face Mask (F3). F3 series have a standard 22 mm tapered that fit external connection to ventilator. This user manual provides you with complete the information you need for the proper use of Full Face Mask (F3).

Intended Use

The Full Face Mask (F3) channels airflow from ventilator systems non-invasively to the patients who suffer from obstructive sleep apnea (OSA) or respiratory insufficiency and respiratory failure.

- The mask is to be used by adult patients (> 30 kg) requiring non-life-support ventilator assistance.
- The mask is intended for single-patient use, for a maximum of seven days, in the hospital or institutional environment.

Before Using the Mask

Full Face Mask (F3) can be matched with four different elbows, leading four different uses situation. See the pictures for details.

CAUTION: It is important to use the correct elbow that matches the respiratory aid system. Please select the correct elbow before using it.

Elbow	V	NV1	NV2	NV3
Picture				
The ventilator systems can be used with	No separate exhalation device required	Require separate exhalation device	Require separate exhalation device	Require separate exhalation device
Color	Transparent	Milky	Blue	Green
Whether have anti-asphyxiation valve (I)	Yes	Yes	No	No
Other feature (s)	/	/	/	Can be used to insert medical pipelines, such as endoscopic ducts (Tracheal intubation with an external diameter less than 8.0 mm)

Contraindications

This mask should be used with caution for patients with following conditions:

- Facial deformity.
- Injured facial skin which will contact with the mask during use.
- Allergic to silicone rubber.

CAUTION: In the US, Federal law restricts the sale of this device to sale by or on the order of a physician.

WARNINGS

- F3 masks are used with elbow "NV1" (main milky), elbow "NV2" (main blue), elbow "NV3" (main green), must have appropriate qualified supervisors, because the patient can not be awake under extreme hypoxia or high carbon dioxide, But unable to take off the mask.

- F3 masks, used with elbow "NV1" (main milky), elbow "NV2" (main blue), elbow "NV3" (main green), can only be used with breathing equipment with alarm function, preventing problems such as low pressure, pipeline shedding, power supply interruption, insufficient gas sources.

- The Exhaust holes (J) of elbow "V" (main transparency) must be kept clear. Occlusion of the exhaust needs to be prevented to avoid having an adverse effect on the safety and quality of the therapy.


- Do not block or seal off the anti-asphyxiation valve when using with elbow "NV1" (main milky) or elbow "NV2" (main blue). Verify that the anti-asphyxiation valve is open when the airflow turned off, the anti-asphyxiation valve is closed when the airflow turned on.

- Use of this mask while the system is not on and operating may cause re-breathing of exhaled air.
- This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.

- To minimize the risk of vomiting during sleep, the patient should avoid eating or drinking for three hours before using the mask. This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Follow all precautions when using supplemental oxygen.

- Oxygen flow must be turned off when the flow generator is not operating, so that unused oxygen does not accumulate within the flow generator enclosure and create a risk of fire.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration varies, depending on the pressure settings, patient breathing pattern, mask, point of application and leak rate.
- The technical specifications of the mask are provided for your clinician to check if it is compatible with the flow generator. If it is used beyond technical specifications or used with incompatible devices, the seal and comfort of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect the function of the flow generator.
- Stop using the Full Face Mask (F3) and consult your physician or sleep therapist, if any of the following symptoms occur: Unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, eye pain, eye infections, tooth, gum, or jaw soreness.
- Refer to your flow generator manual for details on settings and operational information.
- Remove all packaging before using the mask.
- Images shown here are indicative only. If there is inconsistency between the image and actual product, the actual product shall govern.

Fitting the Mask

Note: Taking F3 with elbow "NV1" as an example to show the wearing process.

Fitting the mask tips:

- Check and make sure that the elbow on the mask is the one that the patient needs. Then, pass the upper straps of the headgear through the installation hole of the mask as shown in the figure and secure it, repeat the process with the lower straps. (Fig.1)
- Place the mask on the face and extend the headgear to the back of the head / Please confirm that the label on the back of the headgear faces outward when the headgear is worn. (Fig.2)
- Bring the lower-half straps beneath the ear and encircle the face and then insert it into the fixed hook on the lower-half part of the mask frame. (Fig.3)
- Extend the upper-half straps of the headgear as shown in the above figure until the mask is fixed properly. (Fig.4)
- Repeat the process for the lower-half headgear as shown in the above figure until the mask is fixed properly. (Fig.5)
- Open the device for the air supply. Attach the respiratory hose with the connector on the base of the mask. Adjust the position of the mask to confirm a tight seal. (Fig.6)

CAUTIONS

- The elbow assembly can be dissembled from the frame assembly and changed, but only under the guidance of doctors.
- The elbow assembly of F3 series cannot be disintegrated.
- The valve of elbow "V" and "NV1" cannot be dissembled, or damage may be incurred and the re-installation will become more difficult.

Oxygen / Pressure Port Connector

There are two ports on the base of the mask for taking pressure readings and / or adding additional oxygen.


- Open the Stoppie (F).
- Place the Pressure and / or Oxygen Tube firmly over the port (s).

Note: The two ports are identical; either port may be used for oxygen or pressure measurement. Only remove the Stoppie from ports that are being used.

Cleaning the Mask

This mask is disposable and should be discarded if heavily soiled. If lightly soiled, the mask may be cleaned by wiping with a 70% v/v isopropyl alcohol swab. Do not disassemble any part of the mask except the elbow assembly.

Technical Specifications

Dead Space Information	Dead space is the empty volume of the mask up to the swivel. The dead space of the mask varies according to cushion sizes but is less than 218 mL.
Therapy Pressure	4 to 30 cmH ₂ O
Resistance	Drop in Pressure measured (average for 3 sizes) at 50 L/min: 0.5 cmH ₂ O at 100 L/min: 1.0 cmH ₂ O
Anti-asphyxiation (I) valve closed to Atmosphere (F3 with elbow "V" and "NV1")	≤ 2.8 cmH ₂ O
Anti-asphyxiation (I) valve open to Atmosphere (F3 with elbow "V" and "NV1")	≤ 2.8 cmH ₂ O
Sound (Take the maximum noise of elbow "V")	DECLARED DUAL-NUMBER NOISE EMISSION VALUES in accordance with ISO 4871. The A-weighted sound power level of the mask is no more than 36 dBA. The A-weighted sound pressure level of the mask at a distance of 1 m is no more than 28 dBA.

Notes:

- F3 with elbow "NV2" and "NV3" do not have anti-asphyxiation valve.
- This product is not made with natural rubber latex, PVC or phthalates such as DEHP, DBP or BBP.

Storage

Ensure that the mask is thoroughly clean and dry before storing it for any length of time. Store the mask in a dry place out of direct sunlight.

Disposal

Dispose of in accordance with local regulations.

Symbols

System and Packaging

Caution, consult accompanying documents. Indicates a Warning or Caution and alerts you to a possible injury or explains special measures for the safe and effective use of the device.

Lot number

Catalog number

Not made with natural rubber latex

Temperature limitation

Humidity limitation

Manufacturer

Authorized Representative in the European Community

Single Use only

MANUFACTURER:

BMC MEDICAL CO., LTD.

Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA
 URL: en.bmc-medical.com
 E-mail: intl@bmc-medical.com
 Tel: +86-10-51663880
 Fax: +86-10-51663880 Ext. 810

EU AUTHORISED REPRESENTATIVE:

Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80, 20537 Hamburg, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

Manufacturing Site:

BMC (Tianjin) Medical Co., Ltd.


2/F North Area and 3/F, Building No.4, No.1 Xinxing Road, Wuqing District, (301700) Tianjin, P.R.China
 Tel: +86-22-82939881

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